



- PH & mV Measurement
- 3 1/2 Digit LED display
- Max display 1999 counts
- 0 to 14 PH, 2000 mV
- Temperature compensation
- Two position ABS cabinet

Technical Specifications

DISPLAY	3 1/2 digits LED 0.5" height
RANGE	0 to 14 PH & 0 to 2000 mV DC
RESOLUTION	0.01 pH, 1 mV
ACCURACY	+/- 0.5 % + 5 digits
TEMP. COMPENSATION	0 TO 100 ' C
SELF CHECK	Indicates 7.00 +/- 0.05 pH means the instrument is functioning OK.
ELECTRODE	Provided with the unit, glass combination type
POWER REQUIREMENT	230 +/-10% VAC, 50 Hz, single phase
CABINET SIZE	70 mm (H) x 260 mm (W) x 300 mm (D) Moulded plastic elegant cabinet with two tilt position
WEIGHT	3.0 Kg Approx.
Supplied with glass combination type electrode, Suitable stand, mains cord, spare fuse, and operating manual	

DOT TECHNOLOGIES

(An ISO 9001:2008 Certified Company)

25, Victoria Road, Victoria Layout, Bangalore-560047 INDIA Email: info@dottech.in, web: www.dottech.in

Phone: 080 25549344 / 080 25542753 / 08970010661 / 09844010308 Fax: 080 25542753